# DIAMONDBACK 360° PERIPHERAL PRODUCT REFERENCE GUIDE



© CARDIOVASCULAR SYSTEMS, INC.

CS

# DIAMONDBACK 360° PERIPHERAL ORBITAL ATHERECTOMY SYSTEMS

The CSI<sup>®</sup> Peripheral Orbital Atherectomy System utilizes an eccentrically mounted diamond-coated crown to sand away arterial calcium. As the crown rotates and orbit increases, centrifugal force presses the crown against the lesion, removing or reducing occlusive tissue with each orbit.

#### DIAMONDBACK 360°

PERIPHERAL ORBITAL ATHERECTOMY SYSTEM

MODEL NUMBER	CROWN SIZE	SHAFT LENGTH	QUANTITY	SHEATH SIZE COMPATIBILITY
DIAMONDBACK 360° CLASS	IC CROWN			
DBP-150CLASSIC145	1.50 mm	145 cm	1 each	5 Fr
DBP-200CLASSIC145	2.00 mm	145 cm	1 each	6 Fr
DIAMONDBACK 360° SOLID	CROWN			
DBP-125SOLID145	1.25 mm	145 cm	1 each	4 Fr
DBP-125SOLID200	1.25 mm	200 cm	1 each	5 Fr
DBP-150SOLID145	1.50 mm	145 cm	1 each	5 Fr
DBP-150SOLID200	1.50 mm	200 cm	1 each	5 Fr
DBP-175SOLID180	1.75 mm	180 cm	1 each	5 Fr
DBP-200SOLID145	2.00 mm	145 cm	1 each	6 Fr
DIAMONDBACK 360° MAX C	ROWN			
DBP-200MAX145	2.00 mm	145 cm	1 each	6 Fr
DIAMONDBACK 360° MICRO	CROWN with G	lideAssist®		
DBP-125MICR0145	1.25 mm	145 cm	1 each	4 Fr

### DIAMONDBACK 360°

PERIPHERAL ORBITAL ATHERECTOMY SYSTEM > CEXCHANGEABLE SERIES



SYSTEM MODEL NUMBER Handle + Cartridge	CARTRIDGE MODEL NUMBER Cartridge only	CROWN SIZE	CARTRIDGE SHAFT LENGTH	QUANTITY	SHEATH SIZE COMPATIBILITY
DIAMONDBACK 360° E	Exchangeable Series CLAS	SIC CROW	N with GlideAssis	;†®	
DBP-EX-150CLA145	DBP-CART-150CLA145	1.50 mm	145 cm	1 each	5 Fr
DBP-EX-200CLA145	DBP-CART-200CLA145	2.00 mm	145 cm	1 each	6 Fr
DIAMONDBACK 360° E	Exchangeable Series SOLI		with GlideAssist®		
DBP-EX-125SOL75	DBP-CART-125SOL75	1.25 mm	75 cm	1 each	4 Fr
DBP-EX-125SOL145	DBP-CART-125SOL145	1.25 mm	145 cm	1 each	4 Fr
DBP-EX-125SOL200	DBP-CART-125SOL200	1.25 mm	200 cm	1 each	5 Fr
DBP-EX-150SOL145	DBP-CART-150SOL145	1.50 mm	145 cm	1 each	5 Fr
DBP-EX-150SOL200	DBP-CART-150SOL200	1.50 mm	200 cm	1 each	5 Fr
DBP-EX-175SOL180	DBP-CART-175SOL180	1.75 mm	180 cm	1 each	5 Fr
DBP-EX-200SOL145	DBP-CART-200SOL145	2.00 mm	145 cm	1 each	6 Fr
DIAMONDBACK 360° E	Exchangeable Series MAX	CROWN wit	th GlideAssist®		
DBP-EX-200MAX145	DBP-CART-200MAX145	2.00 mm	145 cm	1 each	6 Fr
DIAMONDBACK 360° E	Exchangeable Series MICR	O CROWN	with GlideAssist®	)	
DBP-EX-125MIC75	DBP-CART-125MIC75	1.25 mm	75 cm	1 each	4 Fr
DBP-EX-125MIC145	DBP-CART-125MIC145	1.25 mm	145 cm	1 each	4 Fr

# ACCESSORY PRODUCTS

Accessory products for use with the CSI Peripheral Orbital Atherectomy Systems to help in the fight against arterial calcium.

#### VIPERWIRE Advance<sup>®</sup> PERIPHERAL GUIDE WIRE WITH FLEX TIP

Multiple enhancements offer improved tip flexibility, navigation, deliverability and trackability compared to our predicate wire. For use with CSI Peripheral Orbital Atherectomy Systems.

	MODEL NUMBER	DESCRIPTION	SHAFT LENGTH	QUANTITY
	VPR-GW-FLEX14	0.014"/0.014" Tip	335 cm	Box of 5
	VPR-GW-FLEX18	0.014"/0.018" Tip	335 cm	Box of 5

VIPERWIRE Advance<sup>®</sup> PERIPHERAL GUIDE WIRE

For use with CSI Peripheral Orbital Atherectomy Systems.

	MODEL NUMBER	DESCRIPTION	SHAFT LENGTH	QUANTITY	
	VPR-GW-14	0.014"/0.014" Tip	335 cm	Box of 5	
uur	VPR-GW-17	0.014"/0.017" Tip	335 cm	Box of 5	-
	VPR-GW-EL14	0.014"/0.014" Tip	475 cm	Box of 5	
	VPR-GW-EL18	0.014"/0.018" Tip	475 cm	Box of 5	

VIPERCATH" XC The longest length exchange catheter available to guide and support guide wire access in the PERIPHERAL EXCHANGE CATHETER vasculature and allow for wire exchanges while providing a conduit for delivering saline and contrast. Designed for kink resistance, torquability and support during delivery. (0.035" guide wire compatible)

MODEL NUMBER	DESCRIPTION	SIZE	LENGTH
035-XC-200-0	Straight	5 Fr	200 cm
035-XC-200-1	30° Angled	5 Fr	200 cm

#### **VIPERTRACK®** RADIOPAQUE TAPE

A radiopaque tape to assist in measuring lesion lengths and marking lesion locations.

1 the	MODEL NUMBER	DESCRIPTION	QUANTITY
Andradan Index 50	VPR-TRK	Radiopaque Tape (30 cm)	Box of 50
	VPR-TRK-20	Radiopaque Tape (30 cm)	Box of 20
	VPR-TRK-MM-20	Radiopaque Tape (500 mm)*	Box of 20

\*With mm markings



An exclusive lubricant designed for use with CSI Orbital Atherectomy Systems.



MODEL NUMBER	DESCRIPTION	QUANTITY
VPR-SLD2	Lubricant	Box of 10

## OAS PUMP

Compact, portable OAS Pump to be used exclusively with CSI Orbital Atherectomy Systems.



MODEL NUMBER	DESCRIPTION	QUANTITY
SIP-3000	OAS Pump	1 System

## **ORDERING INFORMATION**

**Customer Service** 

Hours: 7:00 am - 5:00 pm CST Phone: 877.274.0901 Fax: 612.677.3355 Email: customerservice@csi360.com

The Diamondback 360° Peripheral Orbital Atherectomy System is a percutaneous orbital atherectomy system indicated for use as therapy in patients with occlusive atherosclerotic disease in peripheral arteries and stenotic material from artificial arteriovenous dialysis fistulae. The System is contraindicated for use in coronary arteries, bypass grafts, stents, or where thrombus or dissections are present. Although the incidence of adverse events is rare, potential events that can occur with atherectomy include: pain, hypotension, CVA/TIA, death, dissection, perforation, distal embolization, thrombus formation, hematuria, abrupt or acute vessel closure, or arterial spasm.

Caution: Federal law (USA) restricts this device to sale by, or on the order of, a physician.



1225 Old Hwy 8 NW St. Paul, MN 55112

T: 651.259.1600 877.274.0901

F: 612.677.3355 www.csi360.com

CARDIOVASCULAR SYSTEMS, INC.

For more information, please contact your local CSI representative or call customer service at 1.877.274.0901.

CSI, Diamondback 360, GlideAssist, ViperWire Advance, ViperTrack and ViperSlide are registered trademarks of Cardiovascular Systems. ViperCath is a trademark of Cardiovascular Systems. ©2022 Cardiovascular Systems, Inc. MAT- 0046